

APPLICATION NOTE

[NCMN Central Facility for Scanning Probe and Materials Characterization \(SPMC\)](#)



Electrical Property Mapping at the Nanoscale with PeakForce Tunneling Atomic Force Microscopy (PF-TUNA)

Dr. Lanping Yue, NCMN Specialist

PeakForce TUNA operated in PF-Tapping mode with PF-TUNA Nanoelectrical Application Module provides access to the full fA to μ A current range without changing module or probe holder. PeakForce TUNA enables quantitative conductivity mapping on a wide range of samples while eliminating the adverse effects caused by sample damage and tip contamination. Other benefits include dramatically improved ease of use with ScanAsyst™ Imaging Mode.




Figure 1 High-resolution conductivity map of Co/Pt dots on Chromia. The image size is $5\mu\text{m} \times 5\mu\text{m}$, PF-TUNA images taken at 0.8 V DC bias, using Bruker's SCM-PIT probe (20nm PtIr coating). Image courtesy of Dr. Mahmood, Prof. Binek.




Figure 2

- (a) AFM Topography image
- (b) **TUNA Current** is the average current over one full tapping cycle. It includes both the current measured while tip is in contact with the surface and while it is off the surface.
- (c) **Peak Current** is the instantaneous current, coinciding with Peak Force.

PF-TUNA Remarks

- ✚ Enables reliable nano-electrical imaging on soft delicate samples such as loosely bound nanostructures, conductive polymers
- ✚ Highest resolution current mapping on the most fragile samples such as organic photovoltaics, lithium ion cathodes, and carbon nanotube assemblies
- ✚ Simultaneous nanoelectrical and nanomechanical properties mapping
- ✚ In addition to the imaging mode, PF-TUNA can also measure local current-voltage (I-V) spectra using the spectroscopy mode.